

Jazzy® 623

With features like standard 4-pole motors, Mid-Wheel 6® drive design and ATX Suspension, the Jazzy® 623 is engineered to meet the needs of the most active user. The Jazzy 623 comes standard with high-back seat and Dynamic Shark controller to deliver exceptional performance capability at an excellent value.

*Jazzy® 623
shown in Jazzy Red
with high-back seat*



Pride
Mobility Products
Australia Pty. Ltd.

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



SPECIFICATIONS

Weight Capacity	136 kg (300 lbs.)
Maximum Speed	Up to 8 km/h (5 mph) or up to 9.65 km/h (6 mph)
Ground Clearance	7.62 cm (3") at rear castor
Turning Radius	50.8 cm (20")
Overall Length	90.17 cm (35.5")
Overall Width	61 cm (24")
Drive Wheels	35.56 cm (14") pneumatic or solid
Front Wheels	15.25 cm (6") solid castors
Rear Wheels	20.32 cm (8") solid castors (standard)
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, Mid-Wheel 6®
Braking System	Intelligent braking (electronic, regenerative disc brakes)
Standard Electronics	75 amp, Dynamic Shark controller
Battery Charger	Off-board, 8 amp
Per-Charge Range	Up to 25 km (15.5 miles)
Battery Requirements	(2) 12 volt, deep cycle, 55AH
Battery Weight	17 kg (37 lbs.) each
Base Weight	62 kg (137.5 lbs.)
Standard Seat Weight	20 kg (44 lbs.) (high-back)
Standard Seat	51 x 46-51 cm (20" x 18"-20")
Max. Seat Dimensions	61 x 61 cm (24" x 24")

FEATURES

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® Power Positioning
- Low impact OMNI-Casters (nylon, spherical-shaped castors) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear castor beam for enhanced performance over more varied terrain

OPTIONS & ACCESSORIES*

- Pride® rear bag (includes crutch/ walking stick holder)
- Walker holder (must have Pride® rear bag to attach)
- Oxygen tank holder (must have Pride® rear bag to attach)
- Saddlebag
- Cup holder
- Mobile phone holder

* Please note, accessory bracket may be needed for attaching rear bag.

COLOUR CHOICES



Candy
Apple Red



Viper
Blue

OPTIONAL SEATING

- Contoured high-back seat
- Synergy® Seat
- Specialty Seat
- TRU-Comfort seating
- Quantum dealers only:
TB2 Tilt and Synergy® Seat

Pride[®]
Mobility Products
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

