

Genesis Collection¹



Evolution Collection¹



| MODEL | LC-150S¹ | LC-150M¹ | LC-435S¹ | LC-435M¹ | LC-435L¹ | LC-435XL¹ | LC-435XXL¹ | |
|--------------------------------|------------------------|--|----------------------|--|----------------------|------------------------|---|--|
| Position | Traditional 3-Position | Traditional 3-Position | Enhanced 3-Position | Enhanced 3-Position | Enhanced 3-Position | Traditional 3-Position | Traditional 3-Position | |
| Size | Small | Medium | Small | Medium | Large | Heavy-Duty | Very Heavy-Duty | |
| Suggested Height Range³ | 5'3" & Below | 5'4" - 5'9" | 5'3" & Below | 5'4" - 5'9" | 5'9" - 6' | 5'8" - 6'1" | 5'8" - 6'1" | |
| Weight Capacity | 325 lbs. | 325 lbs. | 400 lbs. | 400 lbs. | 400 lbs. | 500 lbs. | 600 lbs. | |
| Seat to Floor Height | 18" | 20" | 18" | 20" | 22" | 20" | 20" | |
| Seat Depth | 18" | 20" | 18" | 20" | 22" | 22" | 22" | |
| Seat Width | 20" | 22" | 20" | 22" | 24" | 26" | 30" | |
| Top of Back to Seat | 24" | 26" | 24" | 26" | 28" | 29" | 29" | |
| Back Style | Split Back | Split Back | Split Back | Split Back | Split Back | Split Back | Split Back | |
| Fully Padded Chaise | Yes | Yes | Yes | Yes | Yes | Split Ottoman | Split Ottoman | |
| Seating Type | Sinuous Spring | Sinuous Spring | Coil Pocketed Spring | Coil Pocketed Spring | Coil Pocketed Spring | Total Comfort² | Total Comfort² | |
| Distance from Wall | 17" | 18.5" | 16" | 16.5" | 19" | 19" | 19" | |
| Overall Width | 32" | 33" | 38.5" | 40.5" | 42" | 39.5" | 43" | |
| Chair Weight | 111 lbs. | 116 lbs. | 136 lbs. | 141 lbs. | 146 lbs. | 187 lbs. | 196 lbs. | |
| Footrest Extension | N/A | N/A | Standard | Standard | Standard | N/A | N/A | |
| Additional Features | No | No | No | No | No | Lithium Battery Backup | Lithium Battery Backup | |
| Warranty | 2-Yr. Limited | 2-Yr. Limited | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | 7-Yr. Pro-rated | |
| FABRIC SELECTION GUIDE | Standard | Sophie QS LATTE SAPPHIRE SHALE CABERNET | | Cosi QS MOCHA GRAVEL BLUE DIAMOND SYRAH | | | Cosi QS MOCHA GRAVEL | |
| | Optional | Lexis Vinyl QS CHESTNUT | | Sta-Kleen Vinyl QS CHESTNUT | | | Sta-Kleen Vinyl QS CHESTNUT | |

• Recommended HCPCS: E0627—Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.
 • Made of foreign and domestic components



ICON KEY
 = Indicates shortest available lead times

FOOTNOTES:

- Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- Unique spring design and a one-piece soft foam seat for outstanding relaxation.
- Height ranges should be used only as a guide due to varying factors of each individual user.

Like all furniture, Pride® Power Lift Recliners are individually built and there may be differences in actual measurements due to variations in fabric, upholstery, and filling, in addition to the method of measurement.

ADDITIONAL NOTES:

- All chairs feature premium Furniture-Grade Laminate/Hardwood frames.
- All chairs have a Lifetime Warranty on lift mechanism, shell and recline mechanism (except LC-435XXL).
- All chairs have integrated emergency battery backup (batteries not included).
- All chairs can have fiber fill added or removed for custom comfort adjustment.
- All chairs are shipped in an impact-resistant, reinforced, corrugated box.

| POSITION GUIDE | |
|-------------------------------|--|
| PRODUCT VIEW | |
| Traditional 3-Position | The Traditional 3-Position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing. |
| Enhanced 3-Position | The Enhanced 3-Position experience allows the leg rest to extend before the back reclines, offering additional resting positions. |

The lifted, seated and reading positions are achievable on all Pride® Power Lift Recliners.

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

