

MODEL

Jazzy Elite 14-2S-SS ¹	Group 2: Standard (K0822)
Jazzy Elite 14-2S-C ¹	Group 2: Standard (K0823)
Weight Capacity	300 lbs.
Turning Radius	24"
Overall Width ³	24.25"
Overall Length ³	35.5"/42.5" (with footplate)
Maximum Speed ²	Up to 4 mph
Ground Clearance ³	3" at battery tray
Front Wheels	3" solid anti-tips
Drive Wheels	14" foam filled
Rear Wheels	6" solid casters
Drivetrain	Two motor, front-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Limited
High-back Seat Sizes	18" x 18-20", 20" x 18-20" and 22" x 20-22"
Seat-to-Floor Range	20"-22"
Per Charge Range ^{3,5}	Up to 18 miles @ 300 lbs.; Up to 21 miles @ 200 lbs.
Standard Electronics	40A, Dynamic LiNX controller
Component Weights	Base: 99 lbs. (heaviest piece) Seat: 18" wide: 51.8 lbs. 20" wide: 54 lbs. 22" wide: 57.6 lbs. Batteries: 37 lbs. each
Battery Requirements ^{4,6}	(2) 12 volt, NF-22 batteries
Battery Charger	5 amp, off-board
Battery HCPCS Code	E2361
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

⁴ AGM or gel cell type required.

⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁶ Battery weight may vary based on manufacturer.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

0004461-52322

Jazzy[®]
ELITE 14

**EXCELLENT
MANEUVERABILITY
AND OUTDOOR
CAPABILITY**



FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

pridemobility.com |    

TIGHT TURNS AROUND CORNERS

Depth-adjustable,
limited-recline,
fold-flat,
high-back seat

Intuitive controller
is easy to operate

Color-through shroud

Easy rear
access to
batteries

Flip-up foot platform
offers greater comfort

Large 14" drive wheels

FDA Class II
Medical Device¹

The **Jazzy® Elite 14¹** is designed for excellent indoor maneuverability and outstanding outdoor capability. Featuring 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns around corners. Large front wheels offer extra absorption over obstacles for superb climbing capabilities along with front anti-tip wheels for transitioning over grade changes and handling various terrains with ease.

- Color-through shroud keeps the shroud looking new by reducing the visible scrapes and chips
- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel lever
- 60" lap belt

Colors



Jazzy Red

Optional Seating

Deluxe contour, limited recline, high back with solid seat pan, available sizes:

- 16-20"W x 16-20"D

Limited recline high back, available sizes:

- 16"W x 16"D

Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Height adjustable foot platform bracket
- Residual limb support
- Swing-away legrests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger



Front-wheel drive design

Easy-access freewheel levers

40A, Dynamic LiNX controller



Military Patches