

Jazzy

SPORT PORTABLE



The Jazzy® Sport Portable provides a blend of power, performance, style and portability. In-line, front-wheel drive technology provides excellent stability and maneuverability for a solid performance indoors and out while simple disassembly makes it convenient to transport.

CONVENIENT DISASSEMBLY



PERFORMANCE-MATCHED ELECTRONICS



COMFORT HIGH-BACK SEAT WITH HEADREST

AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	25.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Intelligent braking (electronic, regenerative disc brakes)
Standard Electronics	40 amp, PG GC 3 controller
Specialty Controls	N/A
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, deep cycle
Battery Size	U-1
Battery Weight	25 lbs. each
Base Weight	48 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	16" x 16-18", 18" x 18-20", 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy
Apple Red



Viper
Blue

Features

- Easily disassembles into 5 pieces for convenient transport
- Black vinyl, depth adjustable high-back seat
- Flip-up, removable height and width adjustable armrest
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Saddlebag
- XLR USB charger
- Lap belt (50", 70", 80", 90")
- Height-adjustable footrest bracket
- Swing-away footrest
- Elevating leg rests
- Swing-away joystick

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian 260, 440
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder