

Q6 SERIES

THE POWER YOU TRUST



Q6 EDGE®

With features like standard 4-pole motors, Mid-Wheel 6° Drive Design and ATX Suspension, the Q6 Edge® is engineered to meet the performance needs of the most active user. The Q6 Edge accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

Also available with 6 mph high-speed motors

CODE COVERAGE

Q6 EDGE 3S - SS	(Code: K0848)
Q6 EDGE 3S - C	(Code: K0849)
Q6 EDGE 3SP - SS	(Code: K0856)
Q6 EDGE 3SP - C	(Code: K0857)
Q6 EDGE 3MP - SS	(Code: K0861)

Available with
Accu-Trac
Advanced Tracking Technology



Q6000Z

The Quantum® 6000Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6° Design for aggressive outdoor performance and precision indoor maneuverability. The Q6000Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability.

Also available with 8 mph high-speed motors

CODE COVERAGE

Q6000Z 3S - SS	(Code: K0848)
Q6000Z 3S - C	(Code: K0849)
Q6000Z 3SP - SS	(Code: K0856)
Q6000Z 3SP - C	(Code: K0857)
Q6000Z 3MP - SS	(Code: K0861)

Available with
Accu-Trac
Advanced Tracking Technology



Q610

The compact Quantum® 610 is the ultimate choice for users who demand the superb stability of Mid-Wheel 6° Design and tight-quarter maneuverability. Highly adaptable, the Quantum 610 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

CODE COVERAGE


Q610 2S - SS	(Code: K0822)
Q610 2S - C	(Code: K0823)
Q610 2SP - SS	(Code: K0835)
Q610 2SP - C	(Code: K0836)

	Q6 EDGE®	Q6000Z	Q610
Drive Wheels	14" solid or pneumatic ²	14" solid or pneumatic ²	10" solid or pneumatic ²
OMNI-Casters	Front: 6" solid Rear: 6" solid	Front: 6" solid Rear: 6" solid	Front: 5" solid Rear: 6" solid
Suspension	Active-Trac® ATX with enhanced rear suspension	Active-Trac® ATX with integrated rear dual-strut suspension and articulating caster beam	ATX (Active-Trac® with extra stability)
Maximum Speed ¹	Up to 5 mph / Up to 6 mph ²	Up to 6 mph / 8 mph ²	Up to 4 mph
Ground Clearance	3" at center of frame	3.5" at caster arm	2.5" at frame center
Turning Radius	20"	22.5"	21.25"
Overall Size ⁴	Length: 35.5" ¹³ Width: 24"	Length: 36" ¹³ Width: 25.875"	Length: 35.5" ¹³ Width: 23.25"
ALLOWABLE SEAT SIZES	TRU-Balance® 2	TRU-Balance® 3	TRU-Balance® 2
Synergy® Seat	Width: 10" - 22" Depth: 10" - 22"	Width: 12" - 22" Depth: 12" - 22"	Width: 10" - 24" Depth: 10" - 20"
Power Tilt	Width: 14" - 22" Depth: 14" - 22"	Width: 12" - 22" Depth: 12" - 22"	Width: 14" - 20" Depth: 14" - 20"
Lift & Tilt	N.A.	Width: 12" - 22" Depth: 12" - 22"	Width: 16" - 20" Depth: 16" - 20"
Pediatric Tilt	Width: 11" - 14" Depth: 10" - 18"	Width: 11" - 14" Depth: 10" - 18"	Width: 11" - 14" Depth: 10" - 18"
LOWEST SEAT-TO-FLOOR HEIGHTS	TRU-Balance® 3	TRU-Balance® 2	TRU-Balance® 2
Synergy Seat ^{5,6}	15.9" - 17.9"	16.5" - 18.5"	16.5" - 17.5"
Power Tilt ^{5,6}	17.3" - 19.3"	17.5" - 19.5"	16.875"
Lift & Tilt ^{5,6}	N.A.	18" - 20"	18.3"
Pediatric Tilt ^{5,6}	17.25" - 19.25"	17.5" - 19.5"	16.5" - 18.5"
Drivetrain	2-motor, 4-pole, Mid-Wheel 6®	2-motor, Mid-Wheel 6®	2-motor, Mid-Wheel 6®
Batteries	Two 12 volt, deep cycle NF-22	Two 12 volt NF-22 or Group 24	Two 12 volt, deep cycle U-1
Range per Charge ¹	Up to 15.5 miles ¹⁰	Up to 20 miles	Up to 15 miles
Battery Charger	8A, off-board	8A, off-board	5A, off-board
Available Electronics	75A, Q-Logic 2 NE 75A, Q-Logic 2 NE+ 90A, Q-Logic 2 EX	75A, Q-Logic NE 75A, Q-Logic NE+ 90A, Q-Logic EX ¹²	70A PG VR2 80A, Q-Logic EX
Weight Capacity	300 lbs. (user weight)	300 lbs. (user weight)	300 lbs. (user weight) ¹¹
Base Weight	138 lbs.	147 lbs.	112 lbs.
Battery Weight ⁹	38 lbs. (each)	38 lbs. / 52 lbs. (each)	24.5 lbs. (each)
Warranty	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor


FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Optional.
- 3) Requires solid tires.
- 4) Length is measured without front riggings.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Power Recline and Power Front Rigging do not increase seat-to-floor height.
- 7) Up to 24" x 24" with XL versions.
- 8) Required with high-speed motors.
- 9) Battery weight may vary based on manufacturer.
- 10) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4, standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- 11) 250 lbs. weight capacity with TRU-Balance® Power Positioning Systems.
- 12) Required with high speed motors
- 13) Without foot riggings

Q6 Edge® Color Options



Q6000Z Color Options



Q610 Color Options



(US) 866-800-2002 ▪ (Canada) 888-570-1113 ▪ Visit us on the web: www.quantumrehab.com

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.