



EVOLUTION  
LC-435L



**Storage  
Compartment**

*Only available on small,  
medium, and large sizes*




**Cup  
Holder**

*Only available on small,  
medium, and large sizes*

Sitting Pretty™

# Evolution Collection



[PRIDEMOBILITY.COM](http://PRIDEMOBILITY.COM)

Evolution Collection lift chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

While Pride Mobility® began exclusively as a lift chair manufacturer in 1986, our company has never stopped innovating, which makes the **Evolution Collection** a fitting name. This collection fits loungers of all kinds, with sizes ranging from small to XXL.

The small, medium, and large sizes feature a cup holder and storage compartment in the pull-out arms, with USB-A and USB-C chargers in the storage area to keep you connected. These sizes also offer an enhanced three-position experience, which provides additional resting positions.


# Evolution Collection




Model	LC-435S	LC-435M	LC-435L	LC-435XL	LC-435XXL
<b>Size</b>	Small	Medium	Large	X-Large	XX-Large
<b>User Height</b>	5'3" & under	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"
<b>Weight Capacity</b>	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
<b>Seat-to-Floor</b>	18"	20"	22"	20"	20"
<b>Seat Depth</b>	18"	20"	22"	22"	22"
<b>Seat Width</b>	20"	22"	24"	26"	30"
<b>Top of Back to Seat</b>	24"	26"	28"	29"	29"
<b>Fully Padded Chaise</b>	Yes	Yes	Yes	Yes	Yes
<b>Seating Type</b>	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort	Total Comfort
<b>Distance from Wall</b>	16"	16.5"	19"	19"	19"
<b>Overall Width</b>	38.5"	40.5"	42"	39.5"	43"
<b>Total Weight</b>	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
<b>Warranty</b>	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated
<b>Made in America</b> <small>of foreign and domestic parts</small>	No	No	No	Yes	Yes
<b>Position Type</b>	Enhanced 3-position	Enhanced 3-position	Enhanced 3-position	Traditional 3-position	Traditional 3-position

## LC-435S / LC-435M / LC-435L FABRIC OPTIONS

### Cosi Fabrics




Mocha



Gravel



Blue  
Diamond



Syrah

### Sta-Kleen<sup>1</sup>



Chestnut

## FEATURES

- Easy-to-use toggle hand control
- Three-position lift recliners
  - *Enhanced three-position on small, medium, and large for additional resting positions*
- Integrated battery backup in case of power failure
- Chaise pad design fully supports the legs to eliminate pressure points
- Small, medium, and large models feature a cup holder, armrest storage with USB-A/USB-C charging ports, and magnetic hand control dock on both sides
- XL and XXL models feature a USB charger in the hand control, 2-amp lithium battery backup, and are made in the USA of foreign and domestic parts

## OPTIONS & ACCESSORIES

- Lithium battery backup (standard on XL and XXL)
- Head and arm covers
- Non-skid leg levelers
- Cleaning kit

## LC-435XL / LC-435XXL FABRIC OPTIONS

### Cosi Fabrics



Mocha



Gravel

### Sta-Kleen<sup>1</sup>



Chestnut

Evolution Collection lift chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Sitting Pretty™



pridemobility.com

800.800.8586



NOTE: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1) Additional charges will apply for this fabric.