

Pride Lift Chairs

SPECIFICATIONS

	D30	C5	560
Back Width:	74cm (29")	71cm (28")	71cm (28")
Seat To Floor:	48cm (19")	46-51cm (18-20")	48cm (19")
Seat Depth:	56cm (22")	46-51cm (18-20")	48cm (19")
Seat Width:	56cm (22")	56cm (22")	48cm (19")
Top of Back to Seat:	66cm (26")	66-71cm (26-28")	63.5cm (25")
Colour Options:	New Cashmere, New Sapphire, Blush, Forest Green, Levante: Burgundy, Pine Mountain, Royal Blue	Arctic Blue, Birch, Chocolate	Forest Green, Black Cherry, Admiral Blue
	660	805	
Back Width:	71cm (28")	76cm (30")	
Seat To Floor:	48cm (19")	48cm (19")	
Seat Depth:	51cm (20")	51cm (20")	
Seat Width:	51cm (20")	51cm (20")	
Top of Back to Seat:	63.5 (25")	69cm (27")	
Colour Options:	New Cashmere, New Sapphire, Blush	Toffee, Admiral Blue	

BENEFITS

- Engineered furniture laminate/hardwood frame (C5 Steel Frame)
- 170 kg (375 lbs.) weight capacity (136 kg (300 lbs.) weight capacity on TMR805)
- Patented, quiet & smooth lift system
- 2-year limited warranty
- Large side pocket
- Hand control can be easily moved to left hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup

Note: Measurements may vary +/- 1.25 cm (1/2"). Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

38 Lansford Crescent
Avondale, Auckland 1007 New Zealand

Ph: 09 828 0337
Fax: 09 828 0339
www.pridemobility.com

NZ-LC-BIFOLD-10-03-05

Pride
Mobility Products
New Zealand Ltd.



Pride Lift Chairs



The Ultimate in Comfort & Style!

Pride
Mobility Products
New Zealand Ltd.

www.pridemobility.com

Pride Lift Chairs

Pride Lift Chairs

Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any decor, Pride Lift Chairs also feature a 2-year limited warranty. The state-of-the-art lift system utilized in Pride Lift Chairs will recline you further back than any other lift chair, providing wonderful sleeping relaxation. Plus, Pride Lift Chairs are equipped with an emergency battery backup system should the power fail in your home.

DESCRIPTIONS

D30

3-position, medium to large chaise lounger

660

Infinite-position, small to medium chaise lounger with dual motors which enable independent recline/footrest function

C5

3-position, medium to large chaise lounger with fully removable upholstery components

805

2-position, medium to large wall hugger chaise lounger which requires only 10 cm (4") of space between wall and chair

560

3-position, small to medium chaise lounger

PLEASE NOTE: 3-position is a full recline.
2-position is a TV recline.
Infinite-position is a sleep recline.

LIFT CHAIR MODELS



D30

Split-back Design



C5

Sewn Pillow-back Design



560

Pillow-back Design



660

Pillow-back Design




805

Sewn Split-back Design


FABRIC SELECTION ¹




Toffee




Admiral Blue



Black Cherry



Forest




Arctic Blue



Blush




New Cashmere



New Sapphire



Chocolate



Birch

LEVANTE VINYL



Burgundy



Royal Blue



Pine Mountain

¹. Not all fabrics are available for every model. Please see Specifications for details.

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".