

Outstanding Outdoor Performance *with Luxurious Features*



13" Low-profile Tires



Full Suspension



Delta Tiller





Full Lighting Package

The **Pursuit**® incorporates a blend of power, style and performance to handle rugged outdoor terrain. The powerful drivetrain, 24 Volt, 4-pole motor for increased power, 3.56" ground clearance and maximum speed of 8 mph combine power and torque to accelerate your active lifestyle. Standard features include feather touch disassembly, wraparound delta tiller, 13" low-profile solid tires, front and rear suspension, high-back reclining seat, full lighting package, infinitely-adjustable tiller, front basket and rearview mirror. Available in Red and Silver.



Specifications 4-wheel Model # SC713

Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 8 mph
Ground Clearance ²	3.7" at motor, 1.8" at anti-tip, 4.5" at side deck
Turning Radius	67.75" (left) 66.75" (right)
Overall Length ²	53.25"
Overall Width ²	25.25"
Front Tires	3" x 13" solid
Rear Tires	13" solid
Suspension	Front and rear
Range Per Charge ¹	NF-22: Up to 14 miles G24: Up to 19 miles
Total Weight Without Batteries	225 lbs.
Heaviest Piece When Disassembled	95 lbs. (front section)
Standard Seating	Weight: 45.5 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Seat-to-Ground Height	24.75"-26.25"
Seat-to-Deck Height	17"-18.75"
Drivetrain	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery Requirements ^{3,4}	Size: (2) NF-22 (standard) Weight: 38 lbs. each or (2) Group 24 (optional) Weight: 53 lbs. each
Battery Charger	Off-board, 5 amp (NF-22) standard Off-board, 8 amp (Group 24) optional
Colors	  Red Silver

Features

- Easy-to-read LED battery meter
- Reclining high-back seat with sliders, swivel, height adjustable headrest, flip-up, height and width adjustable armrests
- Delta tiller with wraparound handles
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Easy access tiller mounted charge port
- Front and rear suspension
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Full lighting package includes directional signals and large lower headlight for pathway illumination
- Infinitely-adjustable tiller with conveniently located handle
- Large 13" low-profile solid tires
- Powerful drivetrain; features 120 amp microprocessor-based controller
- 24-Volt, 4-pole DC motor for increased power

Accessories

- Cell phone holder
- Cup holder
- 50" Lap belt
- Rear basket
- Rear view mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

Military Patches



LIFT FINDER

Pride offers a variety of lifts to **transport your mobility product with ease**. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC
OR DOES NOT MEET CRITERIA

(US) 800-800-8586

(CAN) 888-570-1113

pridemobility.com

