

3-WHEEL SCOOTER

SPORT RIDER

A STYLISH, COMFORTABLE RIDE

The **Sport Rider** and **Sport Rider Dual Seat** 3-wheel mobility scooters from Pride® combine power, comfort, and style. Both Sport Riders have a complete list of standard features including a fully digital displaying speedometer and odometer, a full suspension system, a high-back adjustable seat (wide bench for Dual Seat), two rearview mirrors, and a powerful front and rear lighting kit. With speeds up to 9.4 mph (up to 7.5 mph for Dual Seat) and a maximum weight capacity of 400 lbs. (Dual Seat accommodates two riders at a weight capacity of 500 lbs.), both Sport Riders have a sleek, sporty look without sacrificing high performance.



Sport Rider with Dual Seat



Windshield
(optional)

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

www.pridemobility.com

SPECIFICATIONS

SPORT RIDER 3-wheel

Weight capacity:	400 lbs., Dual Seat: 500 lbs.
Maximum speed ¹ :	Up to 9.4 mph, Dual Seat up to 7.5 mph
Ground clearance ² :	5"
Turning radius ² :	62.2"
Overall length ² :	65"
Overall width ² :	30.7"
Seat-to-deck height:	18.7"
Seat-to-floor height:	21.8"
Front tires ² :	2.75" w x 14"
Rear tires ² :	2.75" w x 14"
Range per charge ^{1,4} :	Up to 31.25 miles, Dual Seat up to 32.8 miles
Maximum allowable incline:	10°
Total weight without batteries ⁶ :	253.5 lbs., Dual Seat: 257.2 lbs.
Heaviest piece when disassembled:	203.5 lbs. (base weight - does not disassemble)
Standard seating:	Type: Captains, high-back, adjustable, Dual Seat: wide bench Weight: 49.5 lbs., Dual Seat: 53.7 lbs. Material: Black Vinyl Dimensions: 20" x 22", Dual Seat: 36.5" x 17.5"
Motor output:	24V, 1300W, 4300 RPM
Brakes:	Electro-magnetic with drum brake system
Battery requirements ⁵ :	(2) 12-volt, deep cycle Size: (2) 75 AH (standard) or (2) 100AH Weight ³ : (75 AH) 53 lbs. each, (100AH) 75.48 lbs. each
Battery charger:	8-amp, off-board

Color:



Black

FEATURES

- Full digital displaying speedometer/odometer
- Speeds up to 9.4 mph (up to 7.5 for Dual Seat)
- Full suspension system
- 24 volt battery system with (2) 12 volt deep cycle 75 AH batteries (100 AH optional)
- Two rearview mirrors
- Large rear storage pod (optional)
- Windshield (optional)



Full digital displaying speedometer/odometer



Rear storage pod (optional)

1. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

3. Battery weight may vary based on manufacturer.

4. Provided by Wu's Tech. Not tested by Pride.

5. AGM or gel cell type required.

6. Includes standard seat weight.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113
www.pridemobility.com